

Serving
science since
1971



Plastic Containers for Laboratories

FREE Delivery
to mainland
UK



2014



- ***CE marked for in vitro use***
- ***15ml to 16 litres***
- ***Sterile and non-sterile***
- ***With or without labels***
- ***Tamper-evident closures***
- ***From 12p per container***
- ***Jars***
- ***Bottles***
- ***Jerry Cans***
- ***Tubs***
- ***Buckets***
- ***Tubes / Universals***

Contents

Jars

page 4



- polystyrene
- UF or metal screw caps
- 50ml to 300ml
- from 35p

page 6



- polystyrene
- metal or polypropylene screw caps
- 60ml to 250ml
- from 17p

page 8



- polystyrene
- PE screw caps
- 40ml to 200ml
- from 13p

page 10



- polypropylene
- HDPE tamper-evident screw caps
- 25ml to 1000ml
- from 21p

page 12



- polypropylene
- HDPE screw caps
- 150ml to 1000ml
- from 52p

page 14



- polypropylene
- blue or white PE screw caps
- 40ml to 200ml
- from 13p

page 16



- polypropylene
- HDPE tamper-evident screw caps
- 50ml to 1250ml
- from 23p

page 18



- polypropylene
- HDPE screw caps
- 500ml to 1250ml
- from 96p

page 20



- polypropylene
- polypropylene screw caps
- 125ml to 2000ml
- from 54p

page 22



- polypropylene
- HDPE tamper-evident push-on caps
- 350ml to 1000ml
- from 80p

page 24



- PET
- UF or metal screw caps
- 150ml to 300ml
- from 39p

page 26



- PET
- polypropylene screw caps
- 1.6 litres to 4 litres
- from £1.19

Bottles and Jerry Cans

page 28



- PET
- HDPE tamper-evident screw caps
- 250ml to 1000ml. Available dosed
- from 30p

page 30



- PET
- metal or polypropylene tamper-evident screw caps
- 100ml to 1000ml
- from 32p

page 32



- PET or HDPE
- 10 closure types
- 100ml to 1000ml
- from 21p

page 34



- HDPE
- HDPE tamper-evident screw caps
- 150ml to 1000ml
- from 46p

page 36



- HDPE
- HDPE tamper-evident screw caps
- 1 litre to 4 litres
- from £1.07

Tubs and Buckets

page 38



- polypropylene
- tamper-evident snap-on lids
- 250ml and 420ml
- from 25p

page 40



- polypropylene
- polypropylene snap-on lids
- 1 litre to 5 litres
- from 59p

page 42



- polypropylene
- polypropylene tamper-evident snap-on lids
- 1.2 litres to 16.1 litres
- from £1.06

Tubes and Universals

page 44



- polystyrene
- LDPE push-on caps
- 15ml to 165ml
- from 14p

page 46



- polystyrene
- HDPE screw caps
- 30ml. Available bagged
- from 12p

Individually Bagged Jars

page 48



- polystyrene
- metal or UF screw caps
- 60ml to 300ml
- from £1.41

Abbreviations

HDPE	high density polyethylene
LDPE	low density polyethylene
PE	polyethylene
PET	polyethylene terephthalate
UF	urea formaldehyde

About Medfor

Established in 1971, Medfor specialises in the production of plastic sample containers for analytical and clinical use, providing a lighter, more robust and less expensive alternative to glass.

Most products are assembled in-house using a class 10,000 cleanroom, where we cap, label and package containers to provide a wide variety of combinations.

Our products are used by hospitals, government agencies, research companies, food and general industry throughout the UK and other EU countries.

Reliable products

Medfor is registered to ISO 9001:2008 and ISO 13485:2003 and all products in this brochure are CE marked for *in vitro* use.

All our containers are made from food-grade plastics derived from virgin polymers, using components from approved suppliers.

Dependable Service

We provide a friendly and flexible service focussed on meeting your needs.

With a small team of dedicated and experienced staff, you can rely on a familiar voice at the end of the phone and a quick response to emails.

Value for money

Medfor's aim is to offer quality products at competitive prices. Our discount structure provides attractive pricing for both small and large orders, and our minimum order quantity is just one case.

Product Information

Caps

All caps are wadless unless specified.

- **Wadless caps:** Where the cap has been designed to provide a tight seal without the need for a wad.
- **Wadded caps:** A wood pulp or expanded polyethylene (EPE) wad is fitted inside the cap to provide a tighter seal with the neck of the container.

Sterility

Most of our containers are available sterile. This is by gamma irradiation or through being aseptically produced:



Gamma irradiation: This is achieved by a single dose from a radioactive cobalt source (approx 25 - 35 kGy). This is an officially approved process for the sterilisation of single use medical and surgical goods.



Aseptic: Where products are sterilised by the heat at which they are manufactured and then preserved sterile when packaged.

CE marked single use disposable containers

All products are suitable for single *in vitro* use and are CE marked in accordance with the European Directive 98/79/EC.

On labelled containers, the CE mark can be found on the label itself.

For unlabelled products, the CE mark can be found on the outer carton label.

All products are also suitable for general laboratory use.

Easy Ordering

W: www.medfor.co.uk
E: sales@medfor.co.uk
T: +44 (0) 1252 371181
9am to 5pm Monday to Thursday
and 9am to 3pm on Friday
F: +44 (0) 1252 510990

Easy Payment

• We accept the following:



- Proforma: pay by bank transfer. Medfor's bank details will be on your proforma invoice.
- On account: to open a credit account please request an application form.

Free Delivery

Next day delivery is free to all destinations within mainland UK.

Delivery to other destinations is charged at cost. Please check with us before ordering.

Fast Delivery

Many of our products are stock items and will be despatched within 48 hours.

Where it is necessary to assemble an item on receipt of your order, we normally despatch non-sterile products within 3 working days and sterile products within 7 working days.

Stock availability is updated on our website on a daily basis.

Small Orders

For orders of 1 or 2 cases, there is an administration charge of £11.

Orders placed through our website are exempt from this.

Prices

Product prices **include** the cost of the relevant cap, except for the bottles on pages 32-33.

All prices **exclude** VAT. If you are VAT exempt, please include a copy of your exemption certificate with your purchase order.

Prices are in pounds sterling. Whilst we shall endeavour to maintain our prices throughout the life of this brochure, we reserve the right to change them without notice.


Product Specifications

Manufacturers of the containers we supply may change the specifications of the products from time to time. The specifications detailed in this brochure, whether written or pictured, are subject to variation without notice.

Medfor has collated the technical data in this brochure after the exercise of all reasonable care in its compilation, preparation and issue, but it is provided without liability in its application and use.

Clear Polystyrene Honey Jars

- urea formaldehyde screw caps
- metal screw caps



- From 35p
- Ideal for solids, powders, liquids, food and histology samples
- Wide neck facilitates easy filling of liquids or bulky one-piece specimens
- Clarity of plastic provides excellent visibility of contents
- 300ml with metal cap has rubber seal for excellent leak-proof capabilities
- For label designs see page 49
- For individually double-bagged version see page 48

Dimensions


Nominal Capacity (ml)	Brimful Capacity (ml)	Neck Bore (ml)	Overall Diameter (mm)	Overall Height (mm)
50	65	31	43	60
100	115	36	54	70
200	225	54	69	80
300	360	63	80	87

Ordering Information

					Price Per Case			
Code	Capacity (ml)	Sterile ☐	Plain Label	Case Qty	1-4 cases	5-9 cases	10-19 cases	20+ cases
Fitted with wadded white urea formaldehyde screw cap								
HJC0050A	50			315	£146.34	£139.02	£134.63	£128.78
HJC0050B	50		✓	315	£154.70	£146.97	£142.32	£136.14
HJC0050D	50	✓		270	£135.17	£128.41	£124.36	£118.95
HJC0050E	50	✓	✓	270	£142.42	£135.30	£131.03	£125.33
HJC0100A	100			160	£104.37	£99.15	£96.02	£91.85
HJC0100B	100		✓	160	£110.11	£104.60	£101.30	£96.90
HJC0100D	100	✓		150	£106.36	£101.04	£97.85	£93.60
HJC0100E	100	✓	✓	150	£110.39	£104.87	£101.56	£97.14
HJC0200A	200			96	£84.17	£79.96	£77.44	£74.07
HJC0200B	200		✓	96	£86.59	£82.26	£79.66	£76.20
HJC0200D	200	✓		84	£83.13	£78.97	£76.48	£73.15
HJC0200E	200	✓	✓	84	£85.39	£81.12	£78.56	£75.14
HJC0300A	300			60	£62.76	£59.62	£57.74	£55.23
HJC0300B	300		✓	60	£64.35	£61.13	£59.20	£56.63
HJC0300D	300	✓		60	£71.26	£67.70	£65.56	£62.71
HJC0300E	300	✓	✓	60	£72.85	£69.21	£67.02	£64.11
Fitted with wadded white metal screw cap								
HJA0050A	50			315	£124.03	£117.83	£114.11	£109.51
HJA0050B	50		✓	315	£132.40	£125.78	£121.81	£116.51
HJA0050D	50	✓		270	£116.05	£110.25	£106.77	£102.12
HJA0050E	50	✓	✓	270	£123.30	£117.14	£113.44	£108.50
HJA0200A	200			96	£80.44	£76.42	£74.00	£70.79
HJA0200B	200		✓	96	£82.99	£78.84	£76.35	£73.03
HJA0200D	200	✓		84	£79.98	£75.98	£73.58	£70.38
HJA0200E	200	✓	✓	84	£82.24	£78.13	£75.66	£72.37
Fitted with flow-seal metal screw cap								
HJL0300A	300			60	£58.09	£55.19	£53.44	£51.12
HJL0300B	300		✓	60	£59.68	£56.70	£54.91	£52.52
HJL0300D	300	✓		60	£66.59	£63.26	£61.26	£58.60
HJL0300E	300	✓	✓	60	£68.18	£64.77	£62.73	£60.00

Clear Polystyrene Jars

- metal screw caps
- polypropylene screw caps





- From 17p
- Ideal for solids, powders and liquids
- Straight sides facilitate easy and complete recovery of contents
- Clarity of plastic provides excellent visibility of contents
- Metal caps have rubber seals for excellent leak-proof capabilities
- For label designs see page 49
- For individually double-bagged version, see page 48

Dimensions


Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
Fitted with metal screw cap			
60	36	44	61
100	46	55	79
150	46	55	110
250	56	65	119
Fitted with polypropylene screw cap			
60	37	45	61
100	42	51	78

Ordering Information

						Price per case			
Code	Capacity (ml)	Sterile 	Sterile 	Plain Label	Case Qty	1 - 4 cases	5 - 9 cases	10 - 19 cases	20 + cases
Fitted with flow-seal metal screw cap									
SSA0060G	60	✓			250	£53.48	£50.81	£49.20	£47.06
SSA0060H	60	✓		✓	250	£60.43	£57.41	£55.60	£53.18
SSA0060D	60		✓		250	£62.07	£58.97	£57.10	£54.62
SSA0060E	60		✓	✓	250	£69.09	£65.64	£63.56	£60.80
SSA0100G	100	✓			200	£68.45	£65.03	£62.97	£60.24
SSA0100H	100	✓		✓	200	£73.45	£69.78	£67.57	£64.64
SSA0100D	100		✓		175	£68.82	£65.38	£63.31	£60.56
SSA0100E	100		✓	✓	175	£72.62	£68.99	£66.81	£63.91
SSA0150G	150	✓			100	£55.19	£52.43	£50.77	£48.57
SSA0150H	150	✓		✓	100	£58.75	£55.81	£54.05	£51.70
SSA0150D	150		✓		120	£73.83	£70.14	£67.92	£64.97
SSA0150E	150		✓	✓	120	£77.41	£73.54	£71.22	£68.12
SSA0250G	250	✓			50	£37.61	£35.73	£34.60	£33.10
SSA0250H	250	✓		✓	50	£38.86	£36.92	£35.75	£34.20
SSA0250D	250		✓		50	£45.73	£43.44	£42.07	£40.24
SSA0250E	250		✓	✓	50	£47.47	£45.10	£43.67	£41.77
Fitted with white polypropylene screw cap									
SSC0060G	60	✓			50	£49.09	£46.64	£45.16	£43.20
SSC0060H	60	✓		✓	50	£56.03	£53.23	£51.55	£49.31
SSC0060D	60		✓		40	£57.68	£54.80	£53.07	£50.76
SSC0060E	60		✓	✓	40	£64.70	£61.47	£59.52	£56.94
Fitted with wadded white polypropylene screw cap									
SSD0100G	100	✓			30	£76.57	£72.74	£70.44	£67.38
SSD0100H	100	✓		✓	30	£80.29	£76.28	£73.87	£70.66
SSD0100D	100		✓		25	£81.28	£77.22	£74.78	£71.53
SSD0100E	100		✓	✓	25	£85.98	£81.68	£79.10	£75.66

Clear Polystyrene Jars

- polyethylene screw caps




- From 13p
- Ideal for solids, powders and liquids
- Straight sides facilitate easy and complete recovery of contents
- Clarity of plastic provides excellent visibility of contents
- For label design see page 49

Dimensions

Nominal Capacity (ml)	Brimful Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
35	40	28	35	71
55	60	33	40	71
120	125	52	61	74
180	200	52	61	103

Ordering Information

					Price per case			
Code	Capacity (ml)	Sterile 	Plain Label	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Fitted with white polyethylene screw cap								
CPC0040G	35	✓		400	£59.42	£56.45	£54.67	£52.29
CPC0040H	35	✓	✓	400	£75.25	£71.49	£69.23	£66.22
CPC0055G	55	✓		280	£46.82	£44.48	£43.07	£41.20
CPC0055H	55	✓	✓	280	£57.61	£54.73	£53.00	£50.70
CPB0120G	120	✓		380	£155.94	£148.14	£143.46	£137.23
CPB0120H	120	✓	✓	380	£171.22	£162.66	£157.52	£150.67
CPD0200G	180	✓		264	£116.14	£110.33	£106.85	£102.20
CPD0200H	180	✓	✓	264	£126.06	£119.76	£115.98	£110.93

“A wide selection of products and quick despatch.”

Fraser Macgregor, QC Manager
Minelco, Shoreham-by-Sea

Transparent Polypropylene Jars

- HDPE tamper-evident screw caps



- From 21p
- Ideal for solids, powders and liquids
- Wide neck and straight sides facilitate easy and complete recovery of contents
- Transparency of plastic provides good visibility of contents
- Robust plastic minimises likelihood of breakages
- Tamper-evident screw cap provides sample integrity
- For label design see page 49

Dimensions

Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
25	26	32	55
50	35	42	57
100	49	56	63
150	64	72	62
200	64	72	72
250	64	72	87
350	86	95	72
500	86	95	110
1000	108	118	127

Ordering Information

				Price per case			
Code	Capacity (ml)	Plain Label	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Loosely fitted with white HDPE tamper-evident screw cap							
PNC0025A	25		250	£58.96	£56.01	£54.24	£51.88
PNC0025B	25	✓	250	£65.61	£62.33	£60.36	£57.74
PNC0050A	50		250	£74.91	£71.16	£68.92	£65.92
PNC0050B	50	✓	250	£81.55	£77.47	£75.03	£71.76
PNC0100A	100		200	£93.74	£89.05	£86.24	£82.49
PNC0100B	100	✓	200	£99.05	£94.10	£91.13	£87.16
PNC0150A	150		120	£67.56	£64.18	£62.16	£59.45
PNC0150B	150	✓	120	£70.76	£67.22	£65.10	£62.27
PNC0200A	200		96	£65.19	£61.93	£59.97	£57.37
PNC0200B	200	✓	96	£67.74	£64.35	£62.32	£59.61
PNC0250A	250		84	£61.00	£57.95	£56.12	£53.68
PNC0250B	250	✓	84	£63.23	£60.07	£58.17	£55.64
PNC0350A	350		56	£59.65	£56.67	£54.88	£52.49
PNC0350B	350	✓	56	£60.74	£57.70	£55.88	£53.45
PNC0500A	500		36	£43.97	£41.77	£40.45	£38.69
PNC0500B	500	✓	36	£45.10	£42.85	£41.49	£39.69
PNC1000A	1000		18	£43.22	£41.06	£39.76	£38.03
PNC1000B	1000	✓	18	£43.90	£41.71	£40.39	£38.63
As above, in larger case quantities							
PND0350A	350		168	£166.53	£158.20	£153.21	£146.55
PND0350B	350	✓	168	£171.00	£162.45	£157.32	£150.48
PND0500A	500		96	£108.96	£103.51	£100.24	£95.88
PND0500B	500	✓	96	£111.52	£105.94	£102.60	£98.14
PND1000A	1000		60	£130.45	£123.93	£120.01	£114.80
PND1000B	1000	✓	60	£132.05	£125.45	£121.49	£116.20

Transparent Polypropylene Jars

- HDPE screw caps



- From 52p
- Ideal for solids, powders and liquids
- Wide neck and straight sides facilitate easy and complete recovery of contents
- Transparency of plastic provides good visibility of contents
- Robust plastic minimises likelihood of breakages
- For label design see page 49

Dimensions

Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
150	64	72	60
200	64	72	71
250	64	72	86
350	86	95	72
500	86	95	109
1000	108	118	126

Ordering Information


				Price per case			
Code	Capacity (ml)	Plain Label	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Fitted with white HDPE screw cap							
PNA0150A	150		120	£70.58	£67.05	£64.93	£62.11
PNA0150B	150	✓	120	£73.77	£70.08	£67.87	£64.92
PNA0200A	200		96	£67.60	£64.22	£62.19	£59.49
PNA0200B	200	✓	96	£70.15	£66.64	£64.54	£61.73
PNA0250A	250		84	£63.10	£59.95	£58.05	£55.53
PNA0250B	250	✓	84	£65.34	£62.07	£60.11	£57.50
PNA0350A	350		56	£59.65	£56.67	£54.88	£52.49
PNA0350B	350	✓	56	£60.74	£57.70	£55.88	£53.45
PNA0500A	500		36	£43.99	£41.79	£40.47	£38.71
PNA0500B	500	✓	36	£44.87	£42.63	£41.28	£39.49
PNA1000A	1000		18	£46.89	£44.55	£43.14	£41.26
PNA1000B	1000	✓	18	£47.57	£45.19	£43.76	£41.86
As above, in larger case quantities							
PNB0350A	350		168	£166.53	£158.20	£153.21	£146.55
PNB0350B	350	✓	168	£171.00	£162.45	£157.32	£150.48
PNB0500A	500		96	£108.96	£103.51	£100.24	£95.88
PNB0500B	500	✓	96	£111.52	£105.94	£102.60	£98.14
PNB1000A	1000		60	£142.67	£135.54	£131.26	£125.55
PNB1000B	1000	✓	60	£144.27	£137.06	£132.73	£126.96

"Good quality products and the sales and customer team are really efficient."

Max Ferrari, Business Owner
Max 21.Century, Bratislava

Transparent Polypropylene Jars

- Polyethylene screw caps




- From 13p
- Ideal for solids, powders, liquids and food samples
- Distinctive blue version is ideal for use within the food industry
- Wide neck and straight sides facilitate easy and complete recovery of contents
- Robust plastic minimises likelihood of breakages
- Transparency of plastic provides good visibility of contents
- For label design see page 49

Dimensions

Nominal Capacity (ml)	Brimful Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
35	40	28	35	71
55	60	33	40	71
125	130	52	61	74
180	200	52	61	103

Ordering Information

					Price per case			
Code	Capacity (ml)	Sterile 	Plain Label	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Transparent blue jars fitted with blue polyethylene screw cap								
BPC0040G	35	✓		400	£58.91	£55.96	£54.20	£51.84
BPC0040H	35	✓	✓	400	£74.54	£70.81	£68.58	£65.60
BPC0055G	55	✓		280	£44.45	£42.23	£40.89	£39.12
BPC0055H	55	✓	✓	280	£55.09	£52.34	£50.68	£48.48
BPB0125G	125	✓		380	£125.96	£119.66	£115.88	£110.84
BPB0125H	125	✓	✓	380	£141.24	£134.18	£129.94	£124.29
BPD0200G	180	✓		264	£97.87	£92.98	£90.04	£86.13
BPD0200H	180	✓	✓	264	£107.78	£102.39	£99.16	£94.85
Transparent jars fitted with white polyethylene screw cap								
TPC0040G	35	✓		400	£58.95	£56.00	£54.23	£51.88
TPC0040H	35	✓	✓	400	£74.35	£70.63	£68.40	£65.43
TPC0055G	55	✓		280	£42.26	£40.15	£38.88	£37.19
TPC0055H	55	✓	✓	280	£53.05	£50.40	£48.81	£46.68
TPB0125G	125	✓		380	£118.14	£112.23	£108.69	£103.96
TPB0125H	125	✓	✓	380	£133.42	£126.75	£122.75	£117.41
TPD0200G	180	✓		264	£97.07	£92.22	£89.30	£85.42
TPD0200H	180	✓	✓	264	£106.99	£101.64	£98.43	£94.15

*"Very, very helpful, with prompt and efficient service.
Faultless customer service and very friendly."*

Derek Adam, QA Manager
Hyaltech Ltd, Livingston

White Polypropylene Jars

- HDPE tamper-evident screw caps




- From 23p
- Ideal for solids, powders and liquids
- Wide neck and straight sides facilitate easy and complete recovery of contents
- Robust plastic minimises likelihood of breakages
- Tamper-evident screw cap provides sample integrity
- For label design see page 49

Dimensions


Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
50	35	42	57
100	48	56	61
300	72	81	83
500	86	95	108
700	86	95	138
1250	109	118	156

Ordering Information

					Price per case			
Code	Capacity (ml)	Sterile 	Plain Label	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Supplied with white HDPE tamper-evident screw cap – fitted loosely on non-sterile and fitted tightly on sterile products								
WPC0050A	50			250	£66.08	£62.78	£60.79	£58.15
WPC0050B	50		✓	250	£72.80	£69.16	£66.98	£64.06
WPE0050D	50	✓		250	£74.67	£70.94	£68.70	£65.71
WPE0050E	50	✓	✓	250	£81.46	£77.39	£74.94	£71.68
WPC0100A	100			200	£87.95	£83.55	£80.91	£77.40
WPC0100B	100		✓	200	£95.86	£91.07	£88.19	£84.36
WPE0100D	100	✓		200	£96.55	£91.72	£88.83	£84.96
WPE0100E	100	✓	✓	200	£104.51	£99.28	£96.15	£91.97
WPC0300A	300			60	£52.83	£50.19	£48.60	£46.49
WPC0300B	300		✓	60	£54.44	£51.72	£50.08	£47.91
WPE0300D	300	✓		60	£61.43	£58.36	£56.52	£54.06
WPE0300E	300	✓	✓	60	£63.05	£59.90	£58.01	£55.48
WPC0500A	500			36	£42.57	£40.44	£39.16	£37.46
WPC0500B	500		✓	36	£43.70	£41.52	£40.20	£38.46
WPE0500D	500	✓		36	£51.16	£48.60	£47.07	£45.02
WPE0500E	500	✓	✓	36	£52.30	£49.69	£48.12	£46.02
WPC0700A	700			28	£48.26	£45.85	£44.40	£42.47
WPC0700B	700		✓	28	£49.28	£46.82	£45.34	£43.37
WPE0700D	700	✓		28	£56.85	£54.01	£52.30	£50.03
WPE0700E	700	✓	✓	28	£57.88	£54.99	£53.25	£50.93
WPC1250A	1250			15	£38.89	£36.95	£35.78	£34.22
WPC1250B	1250		✓	15	£39.65	£37.67	£36.48	£34.89
WPE1250D	1250	✓		15	£49.70	£47.22	£45.72	£43.74
WPE1250E	1250	✓	✓	15	£48.24	£45.83	£44.38	£42.45
As above, in larger case quantities								
WPD0500A	500			96	£105.23	£99.97	£96.81	£92.60
WPD0500B	500		✓	96	£107.78	£102.39	£99.16	£94.85
WPD0700A	700			80	£128.48	£122.06	£118.20	£113.06
WPD0700B	700		✓	80	£130.69	£124.16	£120.23	£115.01
WPD1250A	1250			45	£105.99	£100.69	£97.51	£93.27
WPD1250B	1250		✓	45	£106.70	£101.37	£98.16	£93.90

White Polypropylene Jars

- HDPE screw caps



- From 96p
- Ideal for solids, powders and liquids
- Wide neck and straight sides facilitate easy and complete recovery of contents
- Robust plastic minimises likelihood of breakages
- For label design see page 49

Dimensions


Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
500	86	95	108
700	86	95	139
1250	109	118	156

Ordering Information

					Price per case			
Code	Capacity (ml)	Sterile ☑	Plain Label	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Fitted with white HDPE screw cap								
WPA0500A	500			36	£42.59	£40.46	£39.18	£37.48
WPA0500B	500		✓	36	£43.47	£41.30	£39.99	£38.25
WPA0500D	500	✓		36	£56.57	£53.74	£52.04	£49.78
WPA0500E	500	✓	✓	36	£56.98	£54.13	£52.42	£50.14
WPA0700A	700			28	£48.26	£45.85	£44.40	£42.47
WPA0700B	700		✓	28	£49.04	£46.59	£45.12	£43.16
WPA0700D	700	✓		28	£56.85	£54.01	£52.30	£50.03
WPA0700E	700	✓	✓	28	£57.65	£54.77	£53.04	£50.73
WPA1250A	1250			15	£41.83	£39.74	£38.48	£36.81
WPA1250B	1250		✓	15	£42.24	£40.13	£38.86	£37.17
WPA1250D	1250	✓		15	£50.42	£47.90	£46.39	£44.37
WPA1250E	1250	✓	✓	15	£50.83	£48.29	£46.76	£44.73
As above, in larger case quantities								
WPB0500A	500			96	£105.23	£99.97	£96.81	£92.60
WPB0500B	500		✓	96	£107.78	£102.39	£99.16	£94.85
WPB0700A	700			80	£128.48	£122.06	£118.20	£113.06
WPB0700B	700		✓	80	£130.46	£123.94	£120.02	£114.80
WPB1250A	1250			45	£114.66	£108.93	£105.49	£100.90
WPB1250B	1250		✓	45	£115.87	£110.08	£106.60	£101.97

White Polypropylene Jars

- polypropylene screw caps



- From 54p
- Ideal for solids and powders
- Wide neck and straight sides facilitate easy and complete recovery of contents
- Robust plastic minimises likelihood of breakages
- For label design see page 49

Dimensions

Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
125	55	61	65
250	64	70	82
350	64	70	115
500	90	98	88
750	108	116	93
1000	108	116	122
2000	130	140	165

Ordering Information

				Price per case			
Code	Capacity (ml)	Plain Label	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Suplied with wadded white polypropylene screw cap							
OWA0125A	125		136	£83.97	£79.77	£77.25	£73.89
OWA0125B	125	✓	136	£87.58	£83.20	£80.57	£77.07
OWA0250A	250		92	£73.22	£69.56	£67.36	£64.43
OWA0250B	250	✓	92	£75.67	£71.89	£69.62	£66.59
OWA0350A	350		60	£58.52	£55.59	£53.84	£51.50
OWA0350B	350	✓	60	£60.13	£57.12	£55.32	£52.91
OWB0500A	500		40	£50.77	£48.23	£46.71	£44.68
OWB0500B	500	✓	40	£51.84	£49.25	£47.69	£45.62
OWB0750A	750		27	£57.90	£55.01	£53.27	£50.95
OWB0750B	750	✓	27	£58.63	£55.70	£53.94	£51.59
OWB1000A	1000		18	£43.68	£41.50	£40.19	£38.44
OWB1000B	1000	✓	18	£44.35	£42.13	£40.80	£39.03
OWA2000A	2000		10	£39.06	£37.11	£35.94	£34.37
OWA2000B	2000	✓	10	£39.48	£37.51	£36.32	£34.74



White Polypropylene Jars

- with polypropylene screw caps
- 70ml, 150ml, 325ml and 415ml capacities
- suitable for autoclaving

For more details, including prices, type "GJA" into the search-box at www.medfor.co.uk

White Polypropylene Jars

- HDPE tamper-evident push-on caps



- From 80p
- Ideal for solids and powders
- Wide neck and straight sides facilitate easy and complete recovery of contents
- Robust plastic minimises likelihood of breakages
- For label design see page 49

Dimensions

Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
350	64	70	115
500	91	96	86
750	109	116	93
1000	109	116	118

Ordering Information

					Price per case			
Code	Capacity (ml)	Sterile ☐	Plain Label	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Supplied with white HDPE tamper-evident push-on cap (fitted on sterile products only)								
WTA0350A	350			60	£54.34	£51.62	£49.99	£47.82
WTA0350B	350		✓	60	£56.17	£53.36	£51.68	£49.43
WTB0350D	350	✓		56	£59.95	£56.95	£55.15	£52.76
WTB0350E	350	✓	✓	56	£61.76	£58.67	£56.82	£54.35
WTA0500A	500			40	£48.95	£46.50	£45.03	£43.08
WTA0500B	500		✓	40	£50.34	£47.82	£46.31	£44.30
WTB0500D	500	✓		36	£53.52	£50.84	£49.24	£47.10
WTB0500E	500	✓	✓	36	£54.87	£52.13	£50.48	£48.29
WTA0750A	750			25	£46.61	£44.28	£42.88	£41.02
WTA0750B	750		✓	25	£47.36	£44.99	£43.57	£41.68
WTB0750D	750	✓		25	£55.20	£52.44	£50.78	£48.58
WTB0750E	750	✓	✓	25	£55.96	£53.16	£51.48	£49.24
WTA1000A	1000			18	£39.75	£37.76	£36.57	£34.98
WTA1000B	1000		✓	18	£40.42	£38.40	£37.19	£35.57
WTB1000D	1000	✓		18	£48.34	£45.92	£44.47	£42.54
WTB1000E	1000	✓	✓	18	£49.02	£46.57	£45.10	£43.14




UN-approved White Polypropylene Jars

- with polypropylene/polyethylene tamper-evident screw caps
- 300ml, 500ml and 1000ml capacities
- air and watertight. UN-X/Y/Z approved for solids and pastes

For more details, including prices, type "CUA" into the search-box at www.medfor.co.uk

Clear PET Jars

- urea formaldehyde screw caps
- metal screw caps




- From 39p
- Ideal for solids, powders, liquids, food and histology samples
- Wide neck facilitates easy filling of liquids or bulky one-piece specimens
- Clarity of plastic provides excellent visibility of contents
- Metal caps have rubber seal for excellent leak-proof capabilities
- For label design see page 49

Dimensions

Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
150	63	76	60
250	63	81	75
300	63	80	87

Ordering Information

					Price per case			
Code	Capacity (ml)	Sterile 	Plain Label	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Fitted with wadded white urea formaldehyde screw cap								
PJB0150A	150			100	£51.65	£49.07	£47.52	£45.45
PJB0150D	150	✓		90	£55.94	£53.14	£51.46	£49.23
PJB0250A	250			60	£35.34	£33.57	£32.51	£31.10
PJB0250B	250		✓	60	£36.94	£35.09	£33.98	£32.51
PJB0250D	250	✓		60	£43.84	£41.65	£40.33	£38.58
PJB0250E	250	✓	✓	60	£45.44	£43.17	£41.80	£39.99
PJB0300A	300			60	£36.34	£34.52	£33.43	£31.98
PJB0300B	300		✓	60	£37.94	£36.04	£34.90	£33.39
PJB0300D	300	✓		60	£44.84	£42.60	£41.25	£39.46
PJB0300E	300	✓	✓	60	£46.44	£44.12	£42.72	£40.87
Fitted with flow-seal metal screw cap								
PJC0150A	150			100	£44.09	£41.89	£40.56	£38.80
PJC0150D	150	✓		90	£48.81	£46.37	£44.91	£42.95
PJC0250A	250			60	£30.67	£29.14	£28.22	£26.99
PJC0250B	250		✓	60	£32.27	£30.66	£29.69	£28.40
PJC0250D	250	✓		60	£39.17	£37.21	£36.04	£34.47
PJC0250E	250	✓	✓	60	£40.76	£38.72	£37.50	£35.87
PJC0300A	300			60	£31.66	£30.08	£29.13	£27.86
PJC0300B	300		✓	60	£33.26	£31.60	£30.60	£29.27
PJC0300D	300	✓		60	£40.16	£38.15	£36.95	£35.34
PJC0300E	300	✓	✓	60	£41.76	£39.67	£38.42	£36.75

Clear PET Jars

- polypropylene screw caps



- From £1.19
- Clarity of plastic provides excellent visibility of contents
- Ideal for larger volumes of liquid samples or storage of food items
- Premium jars with robust walls that are ideal for sample storage and transportation.

Dimensions


Nominal Capacity (L)	Brimful Capacity (L)	Neck Bore (mm)	Overall Width (mm)	Overall Depth (mm)	Overall Height (mm)
1.6	1.8	100	108	108	194
3	3.5	100	150	122	251
4	4.4	100	156	128	292

Ordering Information

			Price per case			
Code	Capacity (L)	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Fitted with white polypropylene screw cap						
SWC1600A	1.6	25	£33.69	£32.01	£30.99	£29.65
SWC3000A	3	20	£30.73	£29.19	£28.27	£27.04
SWC4000A	4	15	£25.60	£24.32	£23.55	£22.53

"From initial order to dispatch of goods and then quick delivery, I am very happy with all aspects of the service provided."

Blaine Wass, Business Owner
BAW Inspection Services, Dronfield



Clear PET Ribbed 5 litre Bottle

- with HDPE tamper-evident screw caps
- neck incorporates carrying handle
- ideal for large volume liquid samples

For more details, including prices, type "RRB" into the search-box at www.medfor.co.uk

Clear PET Square Bottles

- HDPE tamper-evident screw caps




- From 30p
- Ideal for liquid samples, particularly water and beverages
- Clarity of plastic provides excellent visibility of contents
- Square shape minimises storage space
- Lightweight, shatter resistant and durable
- Tamper-evident cap for sample integrity
- Colour coded caps allow for easy identification of samples (White - non-sterile, Green - sterile/non-dosed and blue - sterile/dosed)
- Dosed bottle (1.8% m/v) conforms to Drinking Water Directorate specifications for neutralizing chlorine in water samples
- Dosed versions have the 12 month expiry date clearly printed on each bottle and lot number for complete traceability

Dimensions


Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
250	30	49	142
500	30	63	174
1000	30	74	242

Ordering Information

					Price per case			
Code	Capacity (ml)	Sterile 	Plain Label	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Non-sterile:								
Supplied with white HDPE tamper-evident screw cap								
SQE0250A	250			96	£33.00	£31.35	£30.36	£29.04
SQE0250B	250		✓	96	£35.40	£33.63	£32.57	£31.15
SQE0500A	500			144	£60.32	£57.30	£55.49	£53.08
SQE0500B	500		✓	144	£63.93	£60.73	£58.82	£56.26
SQE1000A	1000			70	£44.13	£41.92	£40.60	£38.83
SQE1000B	1000		✓	70	£45.98	£43.68	£42.30	£40.46
Sterile:								
Fitted with green HDPE tamper-evident screw cap								
SQD0250D	250	✓		84	£37.47	£35.60	£34.47	£32.97
SQD0250E	250	✓	✓	84	£39.75	£37.76	£36.57	£34.98
SQD0500D	500	✓		40	£29.40	£27.93	£27.05	£25.87
SQD0500E	500	✓	✓	40	£30.55	£29.02	£28.11	£26.88
SQD1000D	1000	✓		20	£25.12	£23.86	£23.11	£22.11
SQD1000E	1000	✓	✓	20	£25.66	£24.38	£23.61	£22.58
Dosed:								
Fitted with blue HDPE tamper-evident screw cap and dosed with 1.8% m/v sodium thiosulphate								
SQB0250S	250	✓	dosed	84	£44.21	£42.00	£40.67	£38.90
SQB0500S	500	✓	dosed	40	£32.99	£31.34	£30.35	£29.03
SQB1000S	1000	✓	dosed	20	£27.11	£25.75	£24.94	£23.86

Clear PET Round Bottles

- aluminium screw caps
- polypropylene tamper-evident screw caps




- From 32p
- Ideal for liquid samples, particularly where transportation required
- Robust plastic minimises likelihood of breakages
- Clarity of plastic provides excellent visibility of contents.
- Ordinary screw cap or tamper-evident screw cap which provides sample integrity
- For label design see page 49

Dimensions


Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)	
			Alum. cap	Plastic cap
100	22	48	96	98
250	22	62	129	133
500	22	78	156	161
1000	22	88	219	221

Ordering Information

					Price per case			
Code	Capacity (ml)	Sterile 	Plain Label	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Supplied with white polypropylene tamper-evident screw cap (fitted on sterile products only)								
PTB0100A	100			134	£48.69	£46.26	£44.79	£42.85
PTB0100B	100		✓	134	£52.25	£49.64	£48.07	£45.98
PTC0100D	100	✓		134	£57.19	£54.33	£52.61	£50.33
PTC0100E	100	✓	✓	134	£60.75	£57.71	£55.89	£53.46
PTB0250A	250			60	£31.50	£29.93	£28.98	£27.72
PTB0250B	250		✓	60	£33.33	£31.66	£30.66	£29.33
PTC0250D	250	✓		60	£40.00	£38.00	£36.80	£35.20
PTC0250E	250	✓	✓	60	£41.83	£39.74	£38.48	£36.81
PTB0500A	500			24	£21.59	£20.51	£19.86	£19.00
PTB0500B	500		✓	24	£22.99	£21.84	£21.15	£20.23
PTC0500D	500	✓		24	£30.09	£28.59	£27.68	£26.48
PTC0500E	500	✓	✓	24	£31.28	£29.72	£28.78	£27.53
PTB1000A	1000			15	£18.77	£17.83	£17.27	£16.52
PTB1000B	1000		✓	15	£19.63	£18.65	£18.06	£17.27
PTC1000D	1000	✓		15	£27.27	£25.91	£25.09	£24.00
PTC1000E	1000	✓	✓	15	£28.13	£26.72	£25.88	£24.75
Fitted with wadded aluminium screw cap								
PTA0100A	100			134	£49.50	£47.03	£45.54	£43.56
PTA0100B	100		✓	134	£53.06	£50.41	£48.82	£46.69
PTA0100D	100	✓		134	£57.99	£55.09	£53.35	£51.03
PTA0100E	100	✓	✓	134	£61.56	£58.48	£56.64	£54.17
PTA0250A	250			60	£31.86	£30.27	£29.31	£28.04
PTA0250B	250		✓	60	£33.69	£32.01	£30.99	£29.65
PTA0250D	250	✓		60	£40.36	£38.34	£37.13	£35.52
PTA0250E	250	✓	✓	60	£42.19	£40.08	£38.81	£37.13
PTA0500A	500			24	£21.73	£20.64	£19.99	£19.12
PTA0500B	500		✓	24	£23.09	£21.94	£21.24	£20.32
PTA0500D	500	✓		24	£30.23	£28.72	£27.81	£26.60
PTA0500E	500	✓	✓	24	£31.42	£29.85	£28.91	£27.65
PTA1000A	1000			15	£18.86	£17.92	£17.35	£16.60
PTA1000B	1000		✓	15	£19.72	£18.73	£18.14	£17.35
PTA1000D	1000	✓		15	£27.36	£25.99	£25.17	£24.08
PTA1000E	1000	✓	✓	15	£28.22	£26.81	£25.96	£24.83

Translucent HDPE Bottles / Clear PET Bottles

- choice of closure types



Dimensions

Nominal Capacity (ml)	Brimful Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height Without Closure (mm)
Translucent HDPE Bottles				
250	295	21	56	141
500	540	21	67	186
1000	1100	21	82	243
Clear PET Bottles				
100	120	22	48	95
250	285	22	62	130
500	540	22	79	154
1000	1050	22	88	218



Spray Pump



Spouted Cap



Trigger Spray



Flip Spout



Ribbed Screw Cap



Lock-down Lotion Pump



Child Resistant Closure



Aluminium Screw Cap



Lotion Cap



Tamper-evident Screw Cap

Ordering Information

			Price per case			
Code	Capacity (ml)	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Translucent HDPE bottles only - caps are separate product						
GPA0250A	250	50	£19.52	£18.54	£17.96	£17.18
GPA0500A	500	100	£37.15	£35.29	£34.18	£32.69
GPA1000A	1000	50	£25.89	£24.60	£23.82	£22.78
Clear PET bottles only - caps are separate product						
PTX0100A	100	100	£24.13	£22.92	£22.20	£21.23
PTX0250A	250	50	£21.10	£20.05	£19.41	£18.57
PTX0500A	500	100	£53.12	£50.46	£48.87	£46.75
PTX1000A	1000	50	£40.76	£38.72	£37.50	£35.87

Code	Closure Type	Case Qty	Price per pack*
Packs of closures fit on either type of bottle			
XSP0028A	Spray Pump	50	£17.49
XSC0028A	Spouted Cap	50	£6.57
XTS0028A	Trigger Spray	50	£16.37
XFS0028A	Flip Spout	50	£10.39
XWC0028A	Ribbed Screw Cap	50	£4.78
XLP0028A	Lock-down Lotion Pump	50	£17.01
XCR0028A	Child Resistant Closure	50	£8.93
XAC0028A	Aluminium Screw Cap	50	£6.57
XLC0028A	Lotion Cap	50	£7.75
XTE0028A	Tamper-evident Screw Cap	50	£6.31

- From 21p per bottle
- Ideal for storing and transporting liquids
- Robust plastic minimises likelihood of breakages
- HDPE bottles provide excellent general chemical resistance
- Clarity of PET bottles provides excellent visibility of contents
- Choice of 10 closure types which fit all sizes of either bottle type

***Packs of caps do not contribute to price break quantities**

White HDPE Bottles

- HDPE tamper-evident screw caps




- From 46p
- Ideal for liquid samples, including beverages, chemical and automotive products
- Robust plastic minimises likelihood of breakages
- HDPE provides excellent general chemical resistance
- For label design see page 49

Dimensions

Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
150	25	51	118
250	25	61	131
500	25	77	160
1000	25	97	204

Ordering Information

					Price per case			
Code	Capacity (ml)	Sterile 	Plain Label	Case Qty	1 - 4 cases	5 - 9 cases	10 - 19 cases	20+ cases
Supplied with white HDPE tamper-evident screw cap (fitted on sterile products only)								
STA0150A	150			90	£51.89	£49.30	£47.74	£45.66
STA0150B	150		✓	90	£54.29	£51.58	£49.95	£47.78
STB0150D	150	✓		90	£62.18	£59.07	£57.21	£54.72
STB0150E	150	✓	✓	90	£64.60	£61.37	£59.43	£56.85
STA0250A	250			72	£51.94	£49.34	£47.78	£45.71
STA0250B	250		✓	72	£53.86	£51.17	£49.55	£47.40
STB0250D	250	✓		72	£61.71	£58.62	£56.77	£54.30
STB0250E	250	✓	✓	72	£63.65	£60.47	£58.56	£56.01
STA0500A	500			40	£41.89	£39.80	£38.54	£36.86
STA0500B	500		✓	40	£42.96	£40.81	£39.52	£37.80
STB0500D	500	✓		36	£47.18	£44.82	£43.41	£41.52
STB0500E	500	✓	✓	36	£48.16	£45.75	£44.31	£42.38
STA1000A	1000			18	£27.05	£25.70	£24.89	£23.80
STA1000B	1000		✓	18	£27.54	£26.16	£25.34	£24.24
STB1000D	1000	✓		17	£34.14	£32.43	£31.41	£30.04
STB1000E	1000	✓	✓	17	£34.60	£32.87	£31.83	£30.45
As above, in larger case quantities								
STC0150A	150			300	£155.72	£147.93	£143.26	£137.03
STC0250A	250			225	£146.80	£139.46	£135.06	£129.18
STC0500A	500			100	£95.06	£90.31	£87.46	£83.65
STC1000A	1000			54	£66.95	£63.60	£61.59	£58.92



Translucent HDPE Jars

- with polypropylene screw caps
- 285ml, 525ml, 685ml, 1000ml and 1285ml capacities
- HDPE provides excellent general chemical resistance

For more details, including prices, type "TJA" into the search-box at www.medfor.co.uk

Translucent HDPE Jerry Cans

- HDPE tamper-evident screw caps



- From £1.07
- Ideal for larger volumes of liquid samples e.g. beverage, chemical and automotive products
- Supplied with caps loosely fitted so customer can use tamper-evident feature once jerry can filled
- Robust, shatterproof plastic minimises likelihood of breakages

Dimensions

Nominal Capacity (L)	Brimful Capacity (ml)	Neck Bore (mm)	Overall Width (mm)	Overall Diameter (mm)	Overall Height (mm)
1	1185	31	124	78	155
2	2185	31	137	99	209
4	4300	31	162	127	263

Ordering Information

			Price per case			
Code	Capacity (L)	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Loosely fitted with wadded white HDPE tamper-evident screw cap						
JEA1000A	1	20	£24.34	£23.12	£22.39	£21.42
JEA2000A	2	30	£41.89	£39.80	£38.54	£36.86
JEA4000A	4	20	£33.02	£31.37	£30.38	£29.06

“Medfor provide good quality products that meet our laboratory needs every time we place an order. Reliable service - orders are always processed and despatched quickly.”

Alison Stewart-Gillman
Osmose, Marlow

Clear Polypropylene Round Tubs

- polypropylene tamper-evident snap-on lids



- From 25p
- Ideal for powder, solid and liquid samples, including food, horticulture and paint products
- Stack without lids for efficient storage
- Clarity of plastic provides excellent visibility of contents
- Robust plastic minimises likelihood of breakages

Dimensions

Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
250	83	95	57
420	83	95	93

Ordering Information

			Price per case			
Code	Capacity (ml)	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Supplied with clear polypropylene tamper-evident snap-on lid						
VPA0250A	250	150	£42.05	£39.95	£38.69	£37.00
VPA0420A	420	150	£57.90	£55.01	£53.27	£50.95

"The latest website is uncluttered and easy to use. Packaging is sturdy and the containers inside are always clean and undamaged."

Graham Hodges, Technician
ISP Microcaps (UK) Ltd, Poole



Transparent Polypropylene Cups

- 60ml capacity
- ideal for small volume samples
- wide neck enables easy filling

For more details, including prices, type "S232" into the search-box at www.medfor.co.uk

Translucent Polypropylene Tubs

- polypropylene snap-on lids



- From 59p
- Tubs stack without lids for efficient storage
- Robust plastic minimises likelihood of breakages
- Lids on 1L and 5L are tamper-evident providing sample integrity

Dimensions

Nominal Capacity (L)	Top Inside L x W (mm)	Overall L x W (mm)	Overall Height (mm)
1	176 x 114	185 x 123	72
2	212 x 150	224 x 163	80
3	212 x 150	224 x 163	115
4	212 x 150	224 x 163	149
5	182 x 182	195 x 195	192

Ordering Information

			Price per case			
Code	Capacity (L)	Case Qty	1 - 4 cases	5 - 9 cases	10-19 cases	20+ cases
Supplied with lids						
ICA1000A	1	60	£43.79	£41.60	£40.29	£38.54
ICC2000A	2	20	£27.16	£25.80	£24.99	£23.90
ICC3000A	3	20	£28.90	£27.46	£26.59	£25.43
ICC4000A	4	12	£21.90	£20.81	£20.15	£19.27
ICJ5000A	5	12	£24.76	£23.52	£22.78	£21.79
As above, in larger case quantities						
ICH1000A	1	180	£120.66	£114.63	£111.01	£106.18
ICI2000A	2	80	£93.78	£89.09	£86.28	£82.53
ICI3000A	3	80	£97.99	£93.09	£90.15	£86.23
ICI4000A	4	50	£71.96	£68.36	£66.20	£63.32
ICK5000A	5	48	£81.02	£76.97	£74.54	£71.30

*“Good price and good quality. Strong sturdy pots.
Very knowledgeable team and very helpful every
time.”*

Alex Yem, QA Supervisor
ADM Cocoa UK Ltd, Liverpool

White Polypropylene Buckets

- polypropylene tamper-evident snap-on lids



Dimensions

Brimful Capacity (L)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
1.2	120	139	132
2.7	162	190	157
5.7	208	239	195
10.9	249	283	264
13.2	274	307	257
16.1	311	351	249

Ordering Information

			Price per case			
Code	Capacity (L)	Case Qty	1 - 4 cases	5 - 9 cases	10 -19 cases	20+ cases
Supplied with white polypropylene tamper-evident snap-on lid						
BUA0120A	1.2	10	£12.04	£11.44	£11.08	£10.60
BUA0270A	2.7	10	£16.17	£15.36	£14.88	£14.23
BUA0570A	5.7	10	£19.82	£18.83	£18.23	£17.44
BUA1090A	10.9	10	£26.66	£25.33	£24.53	£23.46
BUA1320A	13.2	10	£39.10	£37.15	£35.97	£34.41
BUA1610A	16.1	10	£43.63	£41.45	£40.14	£38.39


- From £1.06
- Ideal for storing large volumes of liquid
- Handle facilitates easy transportation
- Robust plastic minimises likelihood of breakages and spills

“Friendly, knowledgeable people who remembered me when I phoned to re-order.”

Chris Gray
Mabon Arts, Tredegar

Clear Polystyrene Tubes

- LDPE push-on caps



- From 14p
- Ideal for powders and solids, such as tablets, cake decorations and craft items
- Straight sides permit easy and complete recovery of contents
- Clarity of plastic provides excellent visibility of contents

Dimensions

Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
15	19	23	54
25	25	30	54
45	30	36	70
60	45	50	50
90	45	50	64
130	45	50	92
165	45	50	118

Ordering Information

			Price per case			
Code	Capacity (ml)	Case Qty	1 - 4 cases	5 - 9 cases	10 -19 cases	20+ cases
Supplied with LDPE push-on cap						
FBA0015A	15	200	£31.34	£29.77	£28.83	£27.58
FBA0025A	25	200	£35.55	£33.77	£32.71	£31.28
FBA0045A	45	200	£49.38	£46.91	£45.43	£43.45
FBA0060A	60	100	£33.99	£32.29	£31.27	£29.91
FBA0090A	90	100	£50.12	£47.61	£46.11	£44.11
FBA0130A	130	50	£36.04	£34.24	£33.16	£31.72
FBA0165A	165	50	£42.22	£40.11	£38.84	£37.15

“I find the service to be excellent. The person who handled my order was very courteous and helpful.”

Barbara Watson
Cumbrian Storage Ltd, Workington

Clear Polystyrene Universals

- HDPE screw caps



- From 12p
- For label designs see page 49
- Ideal for patient samples
- Elastomer seal chemically bonded to cap provides excellent seal
- CE marked in accordance with European Directive 98/79/EC
- Intended for *in vitro* use only

Double-bagged version

- Individually double-bagged to preserve the surface sterility of tube until inner bag is opened (inside of tube remains sterile until cap is removed)
- Ideal for hospital theatres and other sterile environments
- Each pack is labelled with expiry date and lot number for full traceability and an irradiation marker to ensure complete sterility

Dimensions

Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
30	22	31	94

Ordering Information

					Price per case			
Code	Sterile A	Sterile R	Label Type	Case Qty	1 - 4 cases	5 - 9 cases	10 -19 cases	20+ cases
Fitted with white HDPE screw cap								
UVA0030G	✓			400	£52.60	£49.97	£48.39	£46.29
UVA0030H	✓		plain	400	£64.05	£60.85	£58.93	£56.36
UVA0030D		✓		400	£63.68	£60.50	£58.59	£56.04
UVA0030E		✓	plain	400	£77.53	£73.65	£71.33	£68.23
Double-bagged with sterility indicator, lot number and expiry date								
UVB0030D-K		✓		80	£123.38	£117.21	£113.51	£108.57
UVB0030F-K		✓	hosp	80	£140.12	£133.11	£128.91	£123.31




Translucent Polypropylene Bijou Tubes

- with polyethylene screw caps
- 7ml capacity
- robust plastic minimises breakages


For more details, including prices, type "BTA" into the search-box at www.medfor.co.uk

Individually Bagged Polystyrene Jars



Dimensions

Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
60	36	44	61
100	46	55	79
150	46	55	110
250	56	65	119




Dimensions

Nominal Capacity (ml)	Neck Bore (mm)	Overall Diameter (mm)	Overall Height (mm)
100	36	54	72
200	54	69	82
300	63	80	89

- Clarity of jars provides excellent visibility of contents
- Wide neck facilitates easy recovery of sample
- Individually double-bagged to preserve surface sterility of jar until inner bag is opened (inside of jar remains sterile until cap is removed)
- Ideal for hospital theatre and other sterile environments
- Each pack is labelled with expiry date and lot number for full traceability and an irradiation marker to ensure complete sterility
- Available with or without 'hospital' label (see opposite for design)
- CE marked in accordance with European Directive 98/79/EC
- Intended for *in vitro* single use only

Ordering Information

					Price per case			
Code	Capacity (ml)	Sterile 	Label Type	Case Qty	1 - 4 cases	5 - 9 cases	10 - 19 cases	20+ cases
Double-bagged with sterility indicator, lot number and expiry date								
SSA0060D-K	60	✓		50	£80.28	£76.27	£73.86	£70.65
SSA0060F-K	60	✓	hosp	50	£90.09	£85.59	£82.88	£79.28
SSA0100D-K	100	✓		40	£69.70	£66.22	£64.12	£61.34
SSA0100F-K	100	✓	hosp	40	£77.75	£73.86	£71.53	£68.42
SSA0150D-K	150	✓		30	£62.51	£59.38	£57.51	£55.01
SSA0150F-K	150	✓	hosp	30	£67.23	£63.87	£61.85	£59.16
SSA0250D-K	250	✓		25	£68.05	£64.65	£62.61	£59.88
SSA0250F-K	250	✓	hosp	25	£63.39	£60.22	£58.32	£55.78
Double-bagged with sterility indicator, lot number and expiry date								
HJC0100D-K	100	✓		50	£103.39	£98.22	£95.12	£90.98
HJC0100F-K	100	✓	hosp	50	£115.99	£110.19	£106.71	£102.07
HJC0200D-K	200	✓		35	£78.69	£74.76	£72.39	£69.25
HJC0200F-K	200	✓	hosp	35	£85.17	£80.91	£78.36	£74.95
HJC0300D-K	300	✓		24	£70.84	£67.30	£65.17	£62.34
HJC0300F-K	300	✓	hosp	24	£72.67	£69.04	£66.86	£63.95


Product Labels (actual size)

Plain labels:

IVD

CE


Ⓢ



IVD

CE

Ⓢ




STERILE

R

IVD

CE

Ⓢ



IVD

CE

Ⓢ

STERILE

A

Hospital label:

IVD

CE


Ⓢ

Name

No./DoB

Ward/GP

Time



STERILE

R

For **free samples** or to place an **order** visit www.medfor.co.uk or call 01252 371181

49

Price Breaks

The price break for an item is based on the aggregate number of cases across a mixed order, not just the quantity for that item.

Eg. For an order of 3 cases of two different items, the '5-9 case' break price will be applied to each, because the aggregate number of cases is six.

Standing Orders

Standing orders can be set up so that pre-arranged quantities of products can be delivered at regular intervals (i.e. weekly, monthly, every 2 months or quarterly).

You benefit from:

- Saving money with better discounts
- Less paperwork since one purchase order covers a whole year's supply
- Peace of mind as the goods are despatched to you automatically
- Being in control with the flexibility to postpone a consignment if you become over-stocked
- Discounts for standing orders are governed by the aggregate number of cases in a consignment, so items in a mixed order will benefit from the overall discount

Please contact us for a quotation or further information.



This brochure contains our most popular products, but there wasn't space for everything in full!

Products with these icons are mentioned in summary format with full details and prices available on our website.



Medfor Products Limited

84a Alexandra Road, Farnborough, GU14 6DD, United Kingdom.

Tel: +44 (0) 1252 371181 Fax: +44 (0) 1252 510990

email: sales@medfor.co.uk www.medfor.co.uk



GB01/52552



GB01/52551